

Drug of abuse saliva test



www.dima-diagnostika.com/oranoxis




The new **Oranoxis**[®] ONSITE Oral-AQ is a quality assured test for drugs of abuse screening. It detects real parent THC and up to 6 additional drugs in saliva within minutes.

Saliva tests have many advantages compared to urine specimen based devices such as making sample adulteration impossible, easy onsite testing anytime anywhere and avoiding urinary retention situations.

Oranoxis[®] saliva tests allow the detection of acute and recent drug usage. With this innovative sample collection system you will collect sufficient saliva for both immediate and confirmatory testing. This extremely simple and user-friendly test is already widely used by:

- Police
- Occupational testing
- Drug substitution centers
- Prisons
- Workplace

Easy to use in 3 steps



ical diagnostics • Oncology • Cardiology • Fertility • Drug of abuse diagnostics • Medical diagnostics • Oncology • Cardiology • Fertility • Drug of abuse diagnostics • Me



DIMA Gesellschaft für Diagnostika mbH • Robert-Bosch-Breite 23 • 37079 Goettingen, Germany
Tel: +49-551-504110 • Fax: +49-551-5041129 • Email: info@dima-diagnostika.com • Website: www.dima-diagnostika.com

Available products

Product name	Description	Content	Product code
Onsite Oral-AQ (7)	AMP / OPI / COC / MET MDMA / BZO / THC	10+1 tests / kit	OR-100704
Onsite Oral-AQ (6)	AMP / OPI / COC / MDMA MET / THC	10+1 tests / kit	OR-100604
Onsite Oral-AQ (THC)	THC	10+1 tests / kit	OR-100100T

Other drug combinations are available on request

Drug detection levels

Drug	Calibrator	Cut-off (ng/ml)
AMP	D-Amphetamine	50
COC	Benzoylcegonine	20
THC	Δ -9-Tetrahydrocannabinol	40
THC	Δ -9-Tetrahydrocannabinol	25
MET	D-Methamphetamine	50
OPI	Opiates	10
BZO	Oxazepam	5
MDMA	\pm 3,4 Methylendioxyamphetamine	75



As a premium manufacturer, Quality and Accuracy are Oranoxis® top priorities:

10+1 quality concept

Oranoxis® kits contain an 11th test for your quality control as well as several quality management documents and a temperature indicator.

Product evaluations

The ONSITE Oral – AQ test has been clinically* evaluated.

Oranoxis® provides high quality innovative products for drugs of abuse, human diseases, food safety and environment testing. We constantly conduct R&D and focus on innovation to develop new products that will best meet customer's needs.

Oranoxis® follows strict quality control practices for our products and our research under Good Manufacturing Practice (GMP), Good Laboratory Practice (GLP) and Quality System Regulation (QSR). Our oral fluid sample collection device is composed of medical grade materials and classified as „Generally Recognized As Safe“ (GRAS) by FDA.

* The study was executed in September 2010 at the Arztpraxis für Medizinische Mikrobiologie (medical practise for medical microbiology), Labordiagnostik & Hygiene Abteilung Drogen-, Medikamenten- und Auftragsanalytik (laboratory diagnostics & hygiene, department drug-, medicine- and analytics on request), Dr. J. Böttcher-Lorenz und Kollegen, Bauhüttenstr.6, 06847 Dessau.